

DemeGUT PLAIN™

Plain Catgut Suture

DemeTech's Plain Catgut Sutures



DemeTech's Plain Catgut suture is an absorbable, sterile, surgical suture composed of highly purified connective tissue (mostly collagen) derived from either beef or sheep intestines. DemeTech's Plain Catgut Sutures are available in Beige. Due to an extremely short absorption period, DemeTech's Plain Catgut elicits only minimal reaction in tissues. The percentage of collagen in the catgut suture often determines the quality of the suture. Higher percentages of collagen allow for: superior tensile strength, longer absorption times, and lower reactions in vivo.

Characteristics of DemeTech's Plain Catgut Sutures:

- DemeTech manufactures its Chromic Catgut sutures with the highest quality collagen, maintaining a composition of 97%-98% of pure collagen.
- DemeTech's Plain Catgut is rapidly absorbed.
- Tensile Strength is maintained for about 10 days.
- Absorption is complete within 70 days.
- DemeTech's Plain Catgut is ideal for tissues which heal rapidly and require minimal support.
- Not CE Approved.

Tensile strength
retained up to 10
days

Tensile Strength

Monofilament



Beige

Material

Complete
Absorption

70 Days

Absorption Profile