

DEMETECH ABSORBABLE SURGICAL SUTURES

SYNTHETIC ABSORBABLE

NATURAL ABSORBABLE

SUTURE STRENGTH RETENTION DAYS



SUTURE BRAND

DemeDIOX™
Polydioxanone

DemeCRYL™
Polyglactin 910

DemeSORB™
Polyglycolic Acid (PGA)

DemeCAPRONE™
Poliglecaprone 25

DemeQUICK™
Rapid Absorbable

DemeGUT™
Chromic Catgut

DemeGUT PLAIN™
Plain Catgut

MATERIAL

Monofilament

Violet Colorless

Multifilament

Violet Colorless

Multifilament

Violet Colorless

Monofilament

Violet Colorless

Multifilament

Violet Colorless

Monofilament

Brown

Monofilament

Beige

TENSILE STRENGTH

14 days: 75% - 80%
28 days: 65% - 70%
42 days: 55% - 60%

7 days: 92%
14 days: 75%
21 days: 50%
28 days: 25%

7 days: 94%
14 days: 84%
21 days: 51%
28 days: 23%

7 days: 70%
14 days: 42%
21 days: 11%

5 days: 50%
10 days: 0%

21-28 days: 0%

7-10 days: 0%

ABSORPTION RATE

182-238 days

56 to 70 days

60 to 90 days

90 to 120 days

42 days

90 days

70 days

BOX



SIZES

7-0 to 2

8-0 to 2

8-0 to 2

7-0 to 2

8-0 to 2

6-0 to 2

6-0 to 2



DemeTECH.
"A Trusted Name in Medicine"

DemeTech Corporation
14175 NW 60th Ave,
Miami Lakes, FL 33014

Phone (305) 824-1048
Fax (305) 437-7607
www.demetech.us