

# DEMETECH ABSORBABLE SURGICAL SUTURES

## SYNTHETIC ABSORBABLE

## NATURAL ABSORBABLE

SUTURE STRENGTH RETENTION DAYS



SUTURE BRAND

**DemeDIOX™**  
Polydioxanone

**DemeCRYL™**  
Polyglactin 910

**DemeSORB™**  
Polyglycolic Acid (PGA)





**DemeCAPRONE™**  
Poliglecaprone 25

**DemeQUICK™**  
Rapid Absorbable

**DemeGUT™**  
Chromic Catgut

**DemeGUT PLAIN™**  
Plain Catgut

MATERIAL



TENSILE STRENGTH

14 days: 75% - 80%  
28 days: 65% - 70%  
42 days: 55% - 60%

7 days: 92%  
14 days: 75%  
21 days: 50%  
28 days: 25%

7 days: 94%  
14 days: 84%  
21 days: 51%  
28 days: 23%

7 days: 70%  
14 days: 42%  
21 days: 11%

5 days: 50%  
10 days: 0%

Tensile strength retained up to 21 days

Tensile strength retained up to 10 days

ABSORPTION RATE

6 Months

56 to 72 days

60 to 90 days

90 to 120 days

42 days

56 to 72 days

70 days

BOX



SIZES

7-0 to 2

8-0 to 2

8-0 to 2

7-0 to 2

8-0 to 2

6-0 to 2

6-0 to 2



**DemeTECH.**  
"A Trusted Name in Medicine"

DemeTECH® Corporation  
14175 NW 60<sup>th</sup> Ave,  
Miami Lakes, FL 33014

Phone: (305) 824-1048  
Fax: (305) 437-7607  
www.demetech.us