

# DEMETECH NON ABSORBABLE SUTURE

## SYNTHETIC SUTURE

## NATURAL SUTURE

**SUTURE BRAND**

**DemeLON™**  
Nylon (Monofilament)

**DemeLON MULTI™**  
Nylon (Multifilament)

**DemeLENE™**  
Polypropylene






**DemeBOND™**  
Polyester (Braided)

**DemeDAC™**  
Polyester (Braided without covering)

**DemeSILK™**  
Silk

**DemeSTEEL™**  
Surgical Steel

**MATERIAL**



**TENSILE STRENGTH**

Gradual loss of tensile strength over time

Gradual loss of tensile strength over time

Indefinite

Indefinite

Indefinite

Gradual loss of tensile strength over time

Indefinite

**BOX**



**SIZES**

10 - 0 to 2

7 - 0 to 2

10 - 0 to 2

8 - 0 to 5

8 - 0 to 5

8 - 0 to 2

8 - 0 to 7



DemeTECH Corporation  
14175 NW 60th Ave,  
Miami Lakes, FL 33014

Phone: (305) 824-1048  
Fax: (305) 437-7607  
www.demetech.us