

DEMETECH NON ABSORBABLE SUTURES

SYNTHETIC SUTURE

NATURAL SUTURE

SUTURE BRAND

DemeLON™
Nylon Monofilament

DemeLON MULTI™
Nylon Multifilament

DemeLENE™
Polypropylene




DemeBOND™
Polyester (Braided)

DemeDAC™
Polyester (Braided
without covering)

DemeSILK™
Silk

DemeSTEEL™
Surgical Steel

MATERIAL



TENSILE STRENGTH

Gradual loss of tensile strength over time

Gradual loss of tensile strength over time

Indefinite

Indefinite

Indefinite

Gradual loss of tensile strength over time

Indefinite

BOX



SIZES

10-0 to 2

7-0 to 2

10-0 to 2

8-0 to 5

8-0 to 5

8-0 to 2

8-0 to 7



DemeTECH®
"A Trusted Name in Medicine"

DemeTECH® Corporation
14175 NW 60th Ave,
Miami Lakes, FL 33014

Phone: (305) 824-1048
Fax: (305) 437-7607
www.demetech.us