

DEMETECH NON ABSORBABLE SUTURE

SYNTHETIC SUTURE

NATURAL SUTURE

SUTURE BRAND

DemeLON™
Nylon Suture (Monofilament)

DemeLON MULTI™
Nylon Suture (Multifilament)

DemeLENE™
Polypropylene Suture

DemeBOND™
Polyester Suture (Braided)

DemeDAC™
polyester Suture (Braided without covering)

DemeSILK™
Silk Suture

DemeSTEEL™
Surgical Steel Suture

MATERIAL

Monofilament

Multifilament

Monofilament

Multifilament

Multifilament

Multifilament

Monofilament

TENSILE STRENGTH

Gradual loss of tensile strength over time

Gradual loss of tensile strength over time

Indefinit

Indefinit

Indefinit

Gradual loss of tensile strength over time

Indefinit

BOX



SIZES

10-0 to 2

7-0 to 2

10-0 to 2

8-0 to 5

8-0 to 5

8-0 to 2

8-0 to 7



DemeTECH
"A Trusted Name in Medicine"

DemeTECH Corporation
14175 NW 60th Ave,
Miami Lakes, FL 33014

Phone: (305) 824-1048
Fax: (305) 437-7607
www.demetech.us