

# DEMETECH NON ABSORBABLE SUTURES

## SYNTHETIC SUTURE

## NATURAL SUTURE

SUTURE BRAND

**DemeLON™**  
Nylon Monofilament

**DemeLON™**  
Nylon Multifilament

**DemeLENE™**  
Polypropylene

**DemeBOND™**  
Polyester (Braided)

**DemeDAC™**  
Polyester (Braided  
without covering)

**DemeSILK™**  
Silk

**DemeSTEEL™**  
Surgical Steel

MATERIAL

Monofilament



Black Clear Blue

Multifilament




Black Blue

Monofilament




Blue

Multifilament



Green White

Multifilament



Green White

Multifilament



Black Blue

Monofilament



Silver

TENSILE STRENGTH

Gradual loss of tensile strength over time

Gradual loss of tensile strength over time

Indefinite

Indefinite

Indefinite

Gradual loss of tensile strength over time

Indefinite

BOX



SIZES

10-0 to 2

7-0 to 2

10-0 to 2

8-0 to 5

8-0 to 5

8-0 to 2

8-0 to 7



**DemeTECH.**  
"A Trusted Name in Medicine"

DemeTech Corporation  
14175 NW 60th Ave,  
Miami Lakes, FL 33014

Phone (305) 824-1048  
Fax (305) 437-7607  
[www.demetech.us](http://www.demetech.us)