


# DemeCAPRONE™

Poliglecaprone 25 Suture


## DemeTech's Poliglecaprone 25 Sutures



DemeTech's fast absorbing Poliglecaprone 25 is made of a synthetic absorbable monofilament composed of glycolide and ε-caprolactone copolymer (PGCL). DemeTech's Poliglecaprone 25 is available in either violet or colorless.


### Characteristics of DemeTech's Poliglecaprone 25 Sutures:

- Very high tensile strength initially
- Excellent handling properties and high pliability
- Very low incidence of infection and trauma
- Smooth tissue passage
- Reliable absorption
- High in vivo strength retention



Tensile Strength

### Monofilament



Material

### Complete Absorption

90 ▶ 120 Days



Absorption Profile