


DemeLON MULTI™

Nylon Suture (Multifilament)

DemeTech's Nylon Sutures



DemeTech's Nylon suture is a non-absorbable, sterile surgical suture composed of the long chain aliphatic polymer Nylon 6.6. DemeTech's Nylon Multifilament sutures is available in black, dyed with FDA listed dye(s). This process is usually followed by gradual encapsulation of the suture by fibrous connective tissue.

Characteristics of DemeTech's Nylon Sutures:

- Composed of tightly braided filaments
- Easy passage through tissues due to its low frictional coefficient.
- Minimal tissue reactivity
- Coated to improve handling.
- Excellent elasticity
- Exceptional for skin closures.

Gradual loss of tensile strength over time

Tensile Strength

Multifilament



Black

Material

Non Absorbable Suture

Absorption Profile