


DemeDAC™

Polyester Suture (Braided without covering)

DemeTech's Polyester Sutures



DemeTech's Polyester (braided without covering) suture is a non-absorbable, sterile, surgical suture composed of Poly (ethylene terephthalate). DemeTech's Polyester Sutures are available in Green and White. DemeTech's Polyester (braided without covering) suture elicits a minimal acute inflammatory reaction in tissues, followed by gradual encapsulation of the suture. High quality materials, combined with precise braiding mechanisms, nearly eliminate the occurrence of postoperative suture fragments in the tissue.

Characteristics of DemeTech's Polyester (braided without covering) Sutures:

- As Polyester (braided without covering) is not coated, it allows for extreme ease of use and secure knot tying.
- It provides for minimal breakage, and thus eliminates the need to remove irritating suture fragments postoperatively (Not CE approved)
- As Polyester (braided without covering) is uncoated, it maintains a higher coefficient of friction as it passes through the tissue.

Indefinite

Tensile Strength

Multifilament



Green

White

Material

Non
Absorbable
Suture

Absorption Profile