

# DemeSTEEL™

Surgical Steel Suture

## DemeTech's Steel Sutures



DemeTech's Surgical Steel sutures are non-absorbable, sterile, surgical sutures composed of high quality stainless steel. Because of its high tensile strength, it exhibits very minimal loss over time and has low tissue reactivity. DemeTech's surgical steel is available as a monofilament.

### Characteristics of DemeTech's Steel Sutures:

- Maximum tensile strength
- Most inert suture

Indefinite

Tensile Strength

Monofilament



Silver

Material

Non  
Absorbable  
Suture

Absorption Profile