

# DemeSTEEL™

## Surgical Steel Suture

### DemeTech's Steel Sutures

DemeTech's Surgical Steel sutures are non-absorbable, sterile, surgical sutures composed of high quality stainless steel. Because of its high tensile strength, it exhibits very minimal loss over time and has low tissue reactivity. DemeTech's surgical steel is available as a monofilament. Surgical Steel is recommended for surgeries ranging from orthopedic surgery to sternal closure. Due to its rigidity, extreme care should be used when handling surgical steel to avoid fragmentation, which could pose harm to both the surgeon and patient.



Indefinite

Tensile Strength

Monofilament



Silver

Material

Non  
Absorbable  
Suture

Absorption Profile

## Characteristics of DemeTech's Steel Sutures:

- Maximum tensile strength
- Most inert suture